



 **ViSafe**™

**REVOLUTIONARY
WEARABLE SENSOR
TECHNOLOGY**


 **dorsaVi**™

DORSAVI

dorsaVi Ltd (ASX: DVL) is an ASX company focused on developing innovative motion analysis device technologies for use in clinical practice, elite sports, and occupational health and safety.

Technology in development since 2000 with commercialization in 2013.

- Patented technology providing accurate movement analysis
- Products across three market segments: Clinical, OH&S and Sport
- FDA Clearance
- Due diligence for clinical market done-validation, reliability, RCT in low back pain




WEARABLE SENSOR TECHNOLOGY



ViSafe Muscle Activity Sensors



ViSafe Movement Sensors



*HOW DOES MOVEMENT
RELATE TO
MUSCULOSKELETAL
INJURIES?*

How does understanding movement help in these 3 sectors:



 **ViMove**™



 **ViPerform**™



 **ViSafe**™

Movement data makes it clear where the focus is needed

GOOD
 MONITOR
 OUTSIDE TEAM NORM

#	Hurdle Step				Overhead Squat		Inline Lunge			
	FX (°)	FX (°)	LFX (°)	LFX (°)	FX (°)	LFX (°)	FX (°)	FX (°)	LFX (°)	LFX (°)
	L	R	L	R	-	-	L	R	L	R
0	11.5	19.1	1.8	-7.4	23.4	-1.8	9	6.3	2.9	3.1
10	16.5	17.8	7.2	-5.8	18.3	0.4	5.8	9.4	3.1	4.9
16	19.8	20.0	7.5	-5.8	21.3	3.5	8.2	8.7	7.1	6.8
49	10.2	12.5	10.4	-6.9	27.3	-0.5	4	3	2.5	1.3
3	8.2	10.0	7.8	-1.9	41.5	-0.2	5.9	4.7	3.5	2.9
5	16.9	20.0	3	-5.6	40.1	-4.9	-4.1	3.3	-4.9	3.9
41	15.3	17.5	10.9	-6.1	19.1	-1	5.5	7.1	5.4	8.2
21	23.8	22.3	10.3	-9.1	31.9	-2	3	5.8	4.1	1.7
50	17.1	14.9	8.2	-9.3	33.8	-1.9	8.7	9.2	1.9	7.4
14	22.2	27.7	9.1	-8.8	19.5	1	7	9	2.3	3.4
15	18.0	18.9	10.9	-5.8	39.9	2.5	6.7	10.5	5	8.6
24	23.2	31.0	6.2	-3.1	28.6	0.9	-5.4	-11.4	-6.4	-3.9
9	18.7	23.2	11.7	-7.1	24.9	3.1	5.3	6.7	1.4	1.4
47	13.0	14.2	11.1	-9.7	15.7	-1.8	3.5	2.6	2.8	4.9
32	17.1	25.6	13.2	-3	23.6	-2.1	-12.5	-6.6	-3.2	-6.7
4	19.1	32.4	6.5	-7.1	34.4	0.4	7.2	11.2	7.6	-6.7
13	19.0	30.4	8.1	-4.9	23.7	-2.2	11.4	8.4	5.4	7.8
34	19.7	13.9	14.1	-11	45	0	8.6	7.7	5.3	7.8
36	21.2	28.3	11.5	-10.5	25.2	-1.2	11.7	15	4.1	3.2
43	15.1	23.4	9.8	-8.5	30.1	-2.8	17.4	11.7	3.7	9.8
11	26.5	32.2	6.3	-7.1	37	2.1	6.2	8.1	4.1	2.6
27	5.1	16.2	7.2	-11.9	32.8	-0.1	9.2	12.6	3.3	10.6
1	14.6	26.0	8	1	34.9	5.5	9.1	17.7	4.2	8.1
48	20.9	23.8	13.9	-4	26	0.6	13.9	10.4	6.8	7.5
44	19.2	21.7	8.2	-9	47.8	-4.3	14.9	9.7	7.2	1.9
37	18.3	25.4	13.5	-8.5	26.3	2.7	6.8	10.1	5	3.9
39	27.6	23.9	13.4	-11.1	22.1	-0.8	6.8	8.9	3.3	10.1
40	28.4	38.5	10.8	-4.3	52.8	0.7	12.6	8.7	4	1.8
38	22.0	20.9	10.1	-7.9	43.4	-1.1	23.3	21	4.8	4.4
23	31.2	24.9	1.1	-2.6	16.1	-1.3	9.9	15.4	6.3	3.5
25	26.0	39.4	7.5	-6.6	42.4	-0.1	12.4	16.9	7	3.2
30	10.6	14.2	15.2	-12.7	29.1	1.5	2.2	8	3.6	10.4
2	19.1	16.1	9.2	-1.4	42.4	1.9	24.1	18.4	2.4	17.1
12	18.4	21.9	15.1	-5.4	32.4	0.6	13.8	22.6	7.2	11.9
17	24.8	24.0	16.5	-0.5	15.5	3.3	19.6	18.4	5.9	8
22	22.2	19.0	15.2	-10.6	52.2	-2	4.5	23.5	4.9	8.6
26	13.0	14.0	14.6	-9.6	47.1	-0.8	14.3	13.5	5.8	14.9
42	28.7	31.0	11.8	-4.7	44.9	-2.6	19.4	11.4	6.9	5.6
46	19.1	27.4	15.4	-4	13.1	2.5	15.9	12.5	5.8	11.2
20	26.5	28.8	15	-9.1	48	1.3	22.2	27	3.9	8.5
31	25.5	39.2	11.7	-4.8	41.7	-0.6	19.1	21.6	4	5.7
33	22.9	30.2	14.6	-1.1	50.2	4.7	16.4	15.2	2.4	5.3
18	30.9	27.5	13.1	-8.6	49	-5	21.2	15.7	-12.9	13.7

Top 3 players played all 22 games



Lower 8 players missed an average of 8 games



HOW DOES MOVEMENT CONTRIBUTE TO WMSDS

Body position

When parts of the body are near the extremes of their range of movements, stretching and compression of tendons and nerves occur. For example sustained forward flexion of spine over 40 degrees


Holding the neck and the shoulders in a fixed position

To perform any controlled movement with the arm, muscles in the shoulder and the neck contract and stay contracted for as long as the task requires.

Vibration

HOW DOES NORMAL & ABNORMAL LOW BACK MOVEMENT LOOK?


Session one



Abnormal Flexion Pattern

- No pelvic movement
- High range of movement

Session eight



Same Patient Regained


- Normal Flexion Pattern
- Normal pelvic movement
- Minimal EMG through flexion range of movement

WHAT DIFFERENCE DOES VISAFE MAKE?


- **Large Casino Group**
 - Reduced injuries on particular role from 72 injuries to 1 after dorsaVi intervention measured the same quarter following year
- **Large retail chain**
 - Since dorsaVi they have experienced a 50% reduction in LTIFRs
- **Large supermarket chain**
 - 87% reduction in manual handling injuries on task assessment post dorsaVi

HOW WE MADE A DIFFERENCE TO THIS CLIENT

- Assessed specific work task
- Used dorsaVi sensors to identify risk & plan changes
- Pre- employment screening
- 300+ screenings performed
- Optimised training protocol – benchmarking



Prior to dorsaVi intervention



After dorsaVi intervention

72 injuries Q1 2014




1 injury Q1 2015

by

- Customised Screening
- Employees suited to task
- Improved work design


MILK CRATE (PROPOSED PROCESS) – 4 BENDS UNLOADING 5 CRATES



MILK CRATE COMPARATIVE ASSESSMENT (SUMMARY)



Woolworths
the fresh food people



- 87% reduction in forward bends > 20° (30 bends versus 4 bends)
- 21% reduction in time (44 secs less time)
- CFO & store manager approval

87% reduction in milk loading injuries in 12 months following trolley change

WHO WE WORK WITH

Manufacturing



Retail



SPECIALTY FASHION | GROUP



Utilities



Transport



Healthcare



Hospitality



Resources



Construction



Others





EVIDENCE COMPARING METHODS OR EQUIPMENT

MINING EXAMPLE

PROPOSED INTERVENTIONS - PRIORITIZATION METHODOLOGY

Impact


- Reduction in manual handling load, force or number of repetitions.
- Improvement in best practice factors.
- Number of personnel or areas affected.
- Savings: costs, time.

Ease of Implementation


- Actions required.
- Costs to implement.
- Time required to implement.

		Ease of implementation	
		Low	High
Impact	Low	LOW	MODERATE
	High	MODERATE	HIGH

PROPOSED INTERVENTIONS - PRIORITIZATION METHODOLOGY



DRIVERS HAVING LOW BACK INJURIES – HOW DO YOU REDUCE RISK?



2A only at a cost of \$7.5M

3D, 4E and 5G at a cost of \$0.5M

THEN HOW DO YOU MAINTAIN LOW INJURY RATES?



Manual Handling Risk Matrix


		Consequence				
		VL - 1	L - 2	M - 3	H - 4	VH - 5
Likelihood	VL - 1	VL 1	VL 2	VL 3	L 4	L 5
	L - 2	VL 2	L 4	L 6	M 8	M 10
	M - 3	VL 3	L 6	M 9	M 12	H 15
	H - 4	L 4	M 8	M 12	H 16	VH 20
	VH - 5	L 5	M 10	H 15	VH 20	VH 25



THEN HOW CAN DORSAVI HELP IMPROVE RTW OUTCOMES

- Workers can wear the sensors for up to 24/7
- Monitoring allows you to analyse what the worker is doing at work or at home that is impacting RTW or CWC
- Biofeedback allows you to program restrictions and alert the worker of 'at risk' movements in real time and real life - guides movement behaviour, avoids aggravation, sustains RTW
- Baseline Assessments can guide when RTW is appropriate
- Monitor and manage treatment effectiveness

DORSAVI TECHNOLOGY ASSISTING RETURN TO WORK





REPORTING

- ViMove delivers easy to read analytics to your desktop in real time
- Assessment Report
 - Normative values
 - Map progression
- Monitoring / Biofeedback

14th August 2013
 record for
Jane Citizen j.citizen7@gmail.com
 practitioner
John Smith

Alert	Possible causes	Suggested management
Poor movement pattern - pelvis > lumbar movement	Reduced lumbar range of motion. maybe easier to strain joint structures due to limited lumbar movement	May benefit from LIVE training/Biofeedback to facilitate single & multi-plane lumbar range of motion, also to monitor change.
Poor movement pattern - pelvis > lumbar movement	Reduced lumbar range of motion. maybe easier to strain joint structures due to limited lumbar movement	May benefit from LIVE training/Biofeedback to facilitate single & multi-plane lumbar range of motion, also to monitor change.





DORSAVI TECHNOLOGY IMPROVING THE HEALTH OF WORKERS

- Objective data for evidence based decision making
- Detailed Task Profiling and Return on Investment Modelling for interventions
- Technology can be used for injury prevention and improving RTW outcomes
- Manual Handling Training that makes a difference
- Experience across a range of industries

